BacterLab Division

Issue date: 02/01/2025

Version: 01.2025

BacterTube™
Bile Esculin Agar (BEA)

BIOCHEMICAL MEDIUM

BacterTube™ Bile Esculin Agar (BEA)

Testing medium used for the isolation and differentiation of *Enterococci*.

Code: 08031

1. INTENDED USE

BacterTube™ Bile Esculin Agar (BEA) is recommended for use as a differential medium in the isolation and presumptive identification of *enterococci*/ group D *streptococci*.

2. PRINCIPLES

BacterTube™ Bile Esculin Agar (BEA) contains esculin, ferric citrate to provide ferric ions, and 4% oxbile to inhibit most other strains of nongroup D *streptococci*. Esculin is hydrolyzed by group D *streptococci* to form dextrose and esculetin. This compound reacts with the ferric ions contained within the medium, turning the medium from its original amber color to a dark brown to black. Thus the tolerance to the presence of oxbile and the hydrolysis of esculin provide the means to presumptively identify group D *streptococci*.

3. TYPICAL COMPOSITION

For 1 liter of medium

Tryptone	17,0 g
Peptic digest of meat	3,0 g
Yeast extract	5,0 g
Bacteriological ox bile	10,0 g
Sodium chloride	5,0 g
Esculin	1,0 g
Ferric ammonium citrate	0,5 g
Sodium azide	0,15 g
Agar	13,0 g

pH of the ready-to-use medium at 25°C: $7,1 \pm 0,2$

4. PREPARATION

The environmental plates are ready-to-use, no preparation required.

5. INSTRUCTIONS FOR USE

- Allow the **BacterTube™ Bile Esculin Agar (BEA)** warm to room temperature before use.
- Inoculate and streak the medium with one isolated pure colony.
- Incubate in an aerobic atmosphere at 35°C for 24-48 hours.
- Observe for growth and blackening of the medium.

NOTE:

- This product is not intended for primary isolation of patient specimens. It should be used only with cultures of isolated organism.
- This product is used in conjunction with other biochemical tests to identify cultures of isolated organism.

6. RESULTS

- The enterococci appear as small translucent colonies surrounded by a black halo.

- Staphylococci and yeasts may yield opaque colonies without a black halo.
- It is indispensable to identify suspected bacteria, especially to eliminate confusion with *Listeria* which may give rise to colonies similar to those of enterococci.

7. QUALITY CONTROL

BacterLab ensures the quality of each product batch by testing with ATCC reference strains.

Reference Strains	Incubation Conditions	Expected Results
<i>Enterococcus faecalis</i> ATCC 29212	35 – 37°C 18 – 24 hours, aerobic	Good growth, black zone around the colony
<i>Streptococcus pyogenes</i> ATCC 19615		Inhibited, no growth

8. STORAGE AND TRANSPORT CONDITIONS

- Storage: 2 – 8°C.
- Transportation: Ambient temperature.

9. PACKAGING

- Packaging: 50 tubes/ box or as per customer request.

10. SHELF LIFE

- Expiration Date: 3 months from the manufacturing date.

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