

BacterLab Division



**Issue date: 02/01/2025**

**Version: 01.2025**

## **MICROBIOLOGICAL CULTURE MEDIUM**

### **BacterPlate™ SDA with Chloramphenicol**

BacterPlate  
SDA có Chloramphenicol

Ready-to-use medium on 90mm plates for selective fungal cultivation. Applied for the selective isolation of fungi from clinical specimens.

**Code: 05028**



## 1. INTENDED USE

**BacterPlate™ SDA with Chloramphenicol** is recommended for the isolation of yeasts and molds, especially when the samples are highly contaminated with bacteria.

The packaging with semi-permeable Cellophane film helps balance the humidity of the environment during storage.

## 2. PRINCIPLES

Peptic digest of Meat is the nitrogen source for growth. Glucose is an energy source.

Chloramphenicol is a heat-stable, broad spectrum antibiotic which inhibits the development of contaminating microflora.

The acid pH favors the growth of yeasts and molds.

## 3. TYPICAL COMPOSITION

*For 1 liter of medium*

Mycological peptone	10,0 g
Glucose	40,0 g
Pancreatic digest of casein	5,0 g
Chloramphenicol	0,008 g
Agar	15,0 g

*pH of the ready-to-use medium at 25°C: 5,6 ± 0,2*

## 4. PREPARATION

The environmental plates are ready-to-use, no preparation required.

## 5. INSTRUCTIONS FOR USE

- Allow the agar plates to reach room temperature. Dry the plates in an incubator by slightly opening the lid.
- Inoculate test samples from the enrichment medium onto the surface of the agar plates.
- Incubate the inoculated plates under aerobic conditions at 25 – 30°C for 72 – 96 hours.

## 6. RESULTS

Separately enumerate yeast colonies and molds

## 7. QUALITY CONTROL

**BacterLab** ensures the quality of each product batch by testing with ATCC reference strains.

Reference Strains	Incubation Conditions	Expected results
<i>C. albicans</i> ATCC 10231	30°C, 3 – 5 days	≥ 50 %
<i>A. brasiliensis</i> ATCC 16404		≥ 50 %
<i>S. cerevisiae</i> ATCC 20827		≥ 50 %
<i>E. coli</i> ATCC 35218		Inhibited

## 8. STORAGE AND TRANSPORT CONDITIONS

- Storage: 2 – 8°C.
- Transportation: Ambient temperature.

## 9. PACKAGING

- Packaging: 10 plates/ box or as per customer request.

## 10. SHELF LIFE

- Expiration Date: 4 months from the manufacturing date.

## 11. BIBLIOGRAPHY

- Solabia Group. *Sabouraud Chloramphenicol Agar*. Biokar Diagnostics. Available at: [https://www.solabia.com/biokar-diagnostics/product/sabouraud-Chloramphenicol-agar/?documentation=1983&\\_wpnonce=66fae9dbeb](https://www.solabia.com/biokar-diagnostics/product/sabouraud-Chloramphenicol-agar/?documentation=1983&_wpnonce=66fae9dbeb).