BacterLab Division

Ngày ban hành: 02/01/2025

Phiên bản: 01.2025

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Polypeptone	23,0 g
Starch	1,0 g
Sodium chloride	5,0 g
Defibrinated sheep blood	50 – 70 mL
Agar	13,5 g

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Polypeptone	23,0 g
Starch	1,0 g
Sodium chloride	5,0 g
Defibrinated horse blood	50 – 70 mL
Bacitracin	1 bộ
Agar	13,5 g

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10. HẠN SỬ DỤNG

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- Neter, E.. 1947. The effect of yeast concentrate on the growth and survival of Haemophilus influenza in infusion broth. Journal of Bacteriology, 54 : 70-71.
- Thayer, J.D. and Martin, J.E.. 1966. Improved medium selective for cultivation of Neisseria gonorrhoeae and Neisseria meningitidis. Public Health Report, 81 : 559-562.
- Illner, P.D., Stoessel, C.J., Drakeford, E. and Vasi, F.. 1966. A new culture medium for medical bacteriology. American Journal of Clinical Pathology, 45 : 502-504.
- NF V08-405. December 1986. Microbiology of food products. Preserves. Detection of {Clostridium} thermophilus..
- NF U47-108. December 2012. Animal health analysis methods - Isolation and identification of Taylorella equigenitalis from Equidae genital sampling.
- NF EN ISO 10272-1. July 2017. Microbiology of the food chain - Horizontal method for detection and enumeration of Campylobacter spp. - Part 1 : Detection method. Modified by NF EN ISO 10272-1/A1 : 2023.
- NF EN ISO 10272-2. July 2017. Microbiology of the food chain — Horizontal method for detection and enumeration of Campylobacter spp. — Part 2: Colony-count technique. Modified by NF EN ISO 10272-2/A1 : 2023.
- NF EN ISO 15213-1. February 2023. Microbiology of the food chain - Horizontal method for the detection and enumeration of Clostridium spp. - Part 1 : enumeration of sulfite-reducing Clostridium spp. by colony-count technique
- European Pharmacopoeia. Chapter 2.6.13. Microbial examination of non-sterile products: test for specified micro organisms. Harmonized method. The United States Pharmacopeia. Chapter Microbiological examination of non-sterile products: Microbial enumeration tests and Chapter Microbiological examination of non-sterile products: Test for specified products.
- The Japanese Pharmacopoeia. Chapter 4.05 Microbial Limit Test I. Microbiological examination of non-sterile products: Total viable aerobic count and II. Microbiological examination of non-sterile products: Test for specified products.
- Solabia Biokar Diagnostics. (2023). Columbia Agar (Base) Technical Data Sheet. Rue des Quarante Mines – ZAC de Ther – Allonne – B.P. 10245 – 60002 Beauvais Cedex, France. Retrieved from <http://www.solabia.com/biokar-diagnostics/>