




SHAKER INCUBATOR KingShak_Incu Elite

- ✓ Simple and efficient operation
- ✓ Convenient and easy to move and clean
- ✓ User-friendly interface
- ✓ 3-year warranty

PURPOSE

- KingShak_Incu Elite is one of the main devices for microbiological culture applications in the medical, pharmaceutical, food, and biological research fields.
- KingShak_Incu Elite is used for shaking under environmental conditions ranging from room temperature + 5°C to 60°C.
- The shaking tray is designed with various options to suit all customer needs.

CONSTRUCTION



SPECIFICATIONS

Model	KingShak_Incu Elite
External dimensions (LxWxH)	L520 x W652 x H485 mm
External dimensions (LxWxH)	L450 x W370 x H26 mm
Oscillation Speeds	5 to 300rpm
Step up the speed	1rpm
Shake amplitude	19mm
Speed stability	±1rpm
Time Range	Up to 99 hour 59 minutes
Maximum Capacity	10kg
Temperature adjustment	from room temperature to 60°C
Temperature adjustment step	1°C
Temperature accuracy	< ±0.5°C
Temperature stability	±0.1°C AT 37°C
Rapid heating time	5°C/minutes
Input voltage	220/240 50hz
Power	55W
Motor	made in Germany
4.3-inch touchscreen display	✓
Made in Vietnam	✓
3-year warranty	✓

ACCESSORIES

Products	Order Code	Description
500mL Erlen Clamp	KingErlenClamp.500	Stainless steel (SUS 304) stand for 500mL Erlen Clamp
1000mL Erlen Clamp	KingErlenClamp.100	Stainless steel (SUS 304) stand for 1000mL Erlen Clamp
96-well Clamp	King96-wellClamp	Stainless steel (SUS 304) stand for 96-well Clamp
Eppendorf Clamp	KingEppendorfClamp	Stainless steel (Inox 304) material can be used for all types of tubes according to customer requirements. The quantity of tubes can be adjusted as needed.
Plate (Optional)		
03 - 500mL Erlen Clamp, 03 - 96-well Clamp	KingPlate.03.03.00	Plate designed to shake 03 - 500mL Erlen Clamp and 03 - 96-well Clamp
03 - 1000mL Erlen Clamp, 02 - 96-well Clamp	KingPlate.03.02.00	Plate designed to shake 03 - 1000mL Erlen Clamp and 02 - 96-well Clamp
06 - 500mL Erlen Clamp	KingPlate.06.00.00	Plate designed to shake 06 - 500mL Erlen Clamp
04 - 1000mL Erlen Clamp	KingPlate.04.00.00	Plate designed to shake 04 - 1000mL Erlen Clamp
06 - 96-well Clamp	KingPlate.00.00.60	Plate designed to shake 06 - 96-well Clamp
Plate for Duran bottle (Optional)		
02 - 1000mL Duran bottle	KingPlateDuran.02.1000	Plate designed to shake 02 - 1000mL Duran bottle
03 - 500mL Duran bottle	KingPlateDuran.03.500	Plate designed to shake 03 - 500mL Duran bottle