



BROCHURE | LabCollect Vacuum Tube



► DESCRIPTION

LabCollect Vacuum Tube® is used as a solution for venous blood collection. LabCollect Vacuum Tube® is used for the transportation and processing of blood samples for serum, plasma, or whole blood testing in clinical laboratories. **LabCollect Vacuum Tube®** are vacuum blood collection tubes with stoppers coded in different colors depending on their function and testing purpose.

► STRUCTURE



► PRODUCTS

LabCollect EDTA K2 Vacuum Tube							
Product Name	Ref.	Volume (mL)	Size (mm)	Description	Material	Label	Cap color
LabCollect EDTA K2 2mL Vacuum tube	1010200	2 mL	13 × 75 mm	EDTAK2	PP	Paper	
LabCollect EDTA K2 3mL Vacuum tube	1010300	3 mL	13 × 75 mm	EDTAK2	PP	Paper	
LabCollect EDTA K2 4mL Vacuum tube	1010400	4 mL	13 × 75 mm	EDTAK2	PP	Paper	
LabCollect EDTA K2 6mL Vacuum tube	1010600	6 mL	13 × 100 mm	EDTAK2	PP	Paper	



LabCollect EDTA K3 vacuum tube							
Product Name	Ref.	Volume (mL)	Size (mm)	Description	Material	Label	Cap color
LabCollect EDTA K3 2mL vacuum tube	1020200	2 mL	13 × 75 mm	EDTAK3	PP	Paper	
LabCollect EDTA K3 3mL vacuum tube	1020300	3 mL	13 × 75 mm	EDTAK3	PP	Paper	
LabCollect EDTA K3 4mL vacuum tube	1020400	4 mL	13 × 75 mm	EDTAK3	PP	Paper	
LabCollect EDTA K3 6mL vacuum tube	1020600	6 mL	13 × 100 mm	EDTAK3	PP	Paper	
LabCollect Lithium Heparin vacuum tube							
Product Name	Ref.	Volume (mL)	Size (mm)	Description	Material	Label	Cap color
LabCollect Lithium Heparin 2mL vacuum tube	1030200	2 mL	13 × 75 mm	Lithium Heparin	PP	Paper	
LabCollect Lithium Heparin 4mL vacuum tube	1030400	4 mL	13 × 75 mm	Lithium Heparin	PP	Paper	
LabCollect Lithium Heparin 6mL vacuum tube	1030600	6 mL	13 × 100 mm	Lithium Heparin	PP	Paper	
LabCollect Sodium Heparin vacuum tube							
Product Name	Ref.	Volume (mL)	Size (mm)	Description	Material	Label	Cap color
LabCollect Sodium Heparin 2mL vacuum tube	1040200	2 mL	13 × 75 mm	Sodium Heparin	PP	Paper	
LabCollect Sodium Heparin 4mL vacuum tube	1040400	4 mL	13 × 75 mm	Sodium Heparin	PP	Paper	
LabCollect Sodium Heparin 6mL vacuum tube	1040600	6 mL	13 × 100 mm	Sodium Heparin	PP	Paper	
LabCollect Tri-Sodium Citrate 3.2% vacuum tube							
Product Name	Ref.	Volume (mL)	Size (mm)	Description	Material	Label	Cap color
LabCollect Tri-Sodium Citrate 3.2% 1,8mL vacuum tube	1060180	1.8 mL	13 × 75 mm	Tri-Sodium Citrate 3.2% 1.8 mL	PET/PP	Paper	
LabCollect Tri-Sodium Citrate 3.2% 2,25mL vacuum tube	1060225	2.25 mL	13 × 75 mm	Tri-Sodium Citrate 3.2% 2.25 mL	PET/PP	Paper	
LabCollect Tri-Sodium Citrate 3.2% 2,7mL vacuum tube	1060270	2.7 mL	13 × 75 mm	Tri-Sodium Citrate 3.2% 2.7 mL	PET/PP	Paper	
LabCollect Tri-Sodium Citrate 3.2% 3,6mL vacuum tube	1060360	3.6 mL	13 × 75 mm	Tri-Sodium Citrate 3.2% 3.6 mL	PET/PP	Paper	
LabCollect Tri-Sodium Citrate 3.2% 4,05mL vacuum tube	1060405	4.05 mL	13 × 75 mm	Tri-Sodium Citrate 3.2% 4.05 mL	PET/PP	Paper	
LabCollect Tri-Sodium Citrate 3.8% vacuum tube							
Product Name	Ref.	Volume (mL)	Size (mm)	Description	Material	Label	Cap color
LabCollect Tri-Sodium Citrate 3.8% 1,8mL vacuum tube	1070180	1.8 mL	13 × 75 mm	Tri-Sodium Citrate 3.8% 1.8 mL	PET/PP	Paper	
LabCollect Tri-Sodium Citrate 3.8% 2,25mL vacuum tube	1070225	2.25 mL	13 × 75 mm	Tri-Sodium Citrate 3.8% 2.25 mL	PET/PP	Paper	
LabCollect Tri-Sodium Citrate 3.8% 2,7mL vacuum tube	1070270	2.7 mL	13 × 75 mm	Tri-Sodium Citrate 3.8% 2.7 mL	PET/PP	Paper	
LabCollect Tri-Sodium Citrate 3.8% 3,6mL vacuum tube	1070360	3.6 mL	13 × 75 mm	Tri-Sodium Citrate 3.8% 3.6 mL	PET/PP	Paper	
LabCollect Tri-Sodium Citrate 3.8% 4,05mL vacuum tube	1070405	4.05 mL	13 × 75 mm	Tri-Sodium Citrate 3.8% 4.05 mL	PET/PP	Paper	
LabCollect Glucose + NaF + K2 Ox vacuum tube							
Product Name	Ref.	Volume (mL)	Size (mm)	Description	Material	Label	Cap color
LabCollect Glucose + NaF + K2 Ox 2mL vacuum tube	1090200	2 mL	13 × 75 mm	NaF + K2 Ox	PET/PP	Paper	
LabCollect Glucose + NaF + K2 Ox 4mL vacuum tube	1090400	4 mL	13 × 75 mm	NaF + K2 Ox	PET/PP	Paper	
LabCollect Glucose + NaF + K2 Ox 5mL vacuum tube	1090500	5 mL	13 × 75 mm	NaF + K2 Ox	PET/PP	Paper	



LabCollect Glucose + NaF + Na Heparin vacuum tube							
Product Name	Ref.	Volume (mL)	Size (mm)	Description	Material	Label	Cap color
LabCollect Glucose + NaF + Na Heparin 2mL vacuum tube	1100200	2 mL	13 × 75 mm	NaF + Na Heparin	PET/PP	Paper	
LabCollect Glucose + NaF + Na Heparin 4mL vacuum tube	1100400	4 mL	13 × 75 mm	NaF + Na Heparin	PET/PP	Paper	
LabCollect Glucose + NaF + EDTA K2 vacuum tube							
Product Name	Ref.	Volume (mL)	Size (mm)	Description	Material	Label	Cap color
LabCollect Glucose + NaF + EDTA K2 2mL vacuum tube	1170200	2 mL	13 × 75 mm	NaF + EDTA K2	PET/PP	Paper	
LabCollect Glucose + NaF + EDTA K2 4mL vacuum tube	1170400	4 mL	13 × 75 mm	NaF + EDTA K2	PET/PP	Paper	
LabCollect Serum vacuum tube							
Product Name	Ref.	Volume (mL)	Size (mm)	Description	Material	Label	Cap color
LabCollect Serum 2mL vacuum tube	1120200	2 mL	13 × 75 mm	Serum	PET/PP	Paper	
LabCollect Serum 3mL vacuum tube	1120300	3 mL	13 × 75 mm	Serum	PET/PP	Paper	
LabCollect Serum 4mL vacuum tube	1120400	4 mL	13 × 75 mm	Serum	PET/PP	Paper	
LabCollect Serum 6mL vacuum tube	1120600	6 mL	13 × 100 mm	Serum	PET/PP	Paper	
LabCollect Serum with gel vacuum tube							
Product Name	Ref.	Volume (mL)	Size (mm)	Description	Material	Label	Cap color
LabCollect Serum with Gel 2.5 mL vacuum tube	1130250	2.5 mL	13 × 75 mm	Serum with gel	PET/PP	Paper	
LabCollect Serum with Gel 3.5 mL vacuum tube	1130350	3.5 mL	13 × 75 mm	Serum with gel	PET/PP	Paper	
LabCollect Serum with Gel 5 mL vacuum tube	1130500	5 mL	13 × 75 mm	Serum with gel	PET/PP	Paper	

► PACKAGING



Packaging Specification
100 tubes/tray
2400 tubes/carton



LABone Scientific Equipment Co.Ltd.,

Address: No. 228/13/3 Nguyen Thi Lang Street, Phu Loi Hamlet, Cu Chi Commune, Ho Chi Minh City, Vietnam
Hotline: 0978 782 147 | Email: info@labone.vn | TAX: 0312088377